


## Avitus<sup>®</sup> Bone Harvester DISTAL FEMUR HARVEST



- 7 minutes harvest time
- 18 CC's of gold-standard cancellous bone
- 10 CC's of gold-standard bone marrow captured
- Graft was successfully used for a Humerus non-union

- A) Avitus<sup>®</sup> cancellous bone retrieved from the Avitus<sup>®</sup> Bone Harvester using the Avitus<sup>®</sup> Filter Insert
- B) Separated Avitus<sup>®</sup> bone marrow retrieved from the Avitus<sup>®</sup> Bone Harvester