



the avitus[®] bone harvester



*The only quick & easy suction
powered debridement
technology*

- Maximize infected bone, osteomyelitis and cancellous bone tumor volume capture.
- Reduce spread of contaminated tissue during debridement.

MINIMIZE SITE ACCESS



I) The Avitus[®] Pilot Hole Creator streamlines cutting an 8 or 11mm cortical window into site.

II) The Avitus[®] Bone Harvester inserts into the cortical hole to begin osteomyelitis debridement.

III) During this process, contaminated contents are collected into the handle of the device.

MAXIMIZE TISSUE CAPTURE

IV) Avitus[®] suction powered technology makes osteomyelitis retrieval quick & easy while maximizing volume of capture. Antibiotic impregnated bone filler can be used to fill defect.